

LATEX-FREE

Dawn

POST-SURGICAL BRA

BY



- 700—Camisole w/2 Removable Drainage Tube Pouches
- 701—w/2 Removable Drainage Tube Pouches and 2 Fiberfill Breast Forms



Dawn




POST-SURGICAL BRA BY WEAR EASE®

The Dawn Bra by Wear Ease® is a post-surgical bra designed for women to wear after breast or chest surgery.


Front View

Inside View

Removable Drainage Tube Pouch



Fiberfill Breast Forms



Inside elastic encircles the entire torso allowing repositioning of the removable drainage tube pouches from front to sides or back according to needs.

The drainage tube pouches can be removed and the camisole can be worn for leisure wear after surgery or during radiation.

AVAILABLE IN SIZES:

Size	XS	S	M	L	XL	1X	2X	3X
Bust	34	36	38	40	42	45	48	51
Waist	26	28	30	32	34	37	40	43
Hip	36	38	40	42	44	47	50	53
Dress size	4/6	8	10	12/14	16/18	20/22	24/26	28

Available in nude and white

DAWN BRA FEATURES:


- Fine, “no-pill” combed cotton
- No inside seams
- Additional under-bust support
- Two Removable Drainage Tube Pouches
- Two triangular fiberfill breast forms (Optional - Style 701)
- Easy to slip-on or step into

Content: 92% combed cotton, 8% spandex jersey.

Washing Instructions: Machine wash gentle cold, use non-chlorine bleach if needed. Do not use bleach or fabric softener. Tumble dry low or hang to dry.

Remove the breast forms before washing. When washing the breast forms, remove the fiberfill. The pouches can be washed attached to the camisole or washed separately.

AVAILABLE IN COLORS:



Phone: 866.251.0076
Fax: 208.445.0957
www.wearease.com

